

FM-02-124-C

2023 2 2  
2015 3 1

GB/T 16483 GB/T 17519  
19375  
4.0

---

o-Phthalic anhydride

295

515000

jhd-service@ghitech.com

0754-88221999

0754-82515813

0532-83889090



**GHS**

, ( 4)  
( 1)  
( 1)  
( 1)  
( 1)  
( 3)



295

**E-mail** jhd-service@ghitech.com

1 / 7

(86) 0754-88213888

**www.ghitech.com**

●

P260 / / / / /

P264

P270

P271

P272

P280 / /

P285

●

P301 + P312

P302 + P352

P304 + P340 ,

P305 + P351 + P338

P310

P330

P333 + P313 /

P362

●

P405

●

P501 /

/

%

99.7

**CAS**

85-44-9



129~133  
285  
152  
[%  
( =1) 1.53



**PH**

=1 5.10

**kJ/mol**

[% ] 1.7

**kPa**

**MPa**

/

IARC:

0.1%

IARC

-

(LC50) - - > 99 mg/l - 96 h

EC50 - Daphnia magna ( ) - 71 mg/l - 48 h

ErC50 - Pseudokirchneriella subcapitata - 68 mg/l - 72 h

EC50 - - > 1,000 mg/l - 3 h

/ - 14 d

: 85 % -

UN 2214

8

III

-

/

-

FE

